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All products manufactured by Felisa are guaranteed for one year from date of shipment from Felisa plant. All those products returned within one year will be rebuild or replaced under the guarantee regardless of reason for failure. Electric motors, are guaranteed according to the manufacturer policies.

Exceptions may be made by Felisa on particular applications, where experience has indicated conditions are so unusual that premature failure can be expected. Transportations on charges, in all cases will be at customer expense. Maximum liability is in no case to exceed the value of the Felisa product involved.

Felisa has the right to change or modify the products in order to improve their presentation and /or operation.

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
(14)

REPAIRS.

Like any other manufactured product, some parts of the Oven could be damaged after using for a long time. To replace them always use original factory parts, a list of which is included in this instruction manual. Always order the parts with its corresponding number. All the parts can be ordered from any distributor or directly from Felisa.

HEATING ELEMENT.

Nickel-Chrome heating element on stainless steel tubing. This type of heating element supplies a high efficiency in heat radiation.

It should be taken the precaution of no causing damage to the element with any tool inserted when operating the Oven or spilling liquids inside the chamber. All the heating elements should be considered perishable and therefore replaceable, however, a reasonable care in its use, will greatly extend the service life. Since the manufacturer does not have control over the use and care of these elements, no guarantee can be made.

IMPORTANT.

To carry out any maintenance work, disconnect the Incubator from its source of energy.

Do not Spill any solution inside the camber.

Do not change the position of the temperature sensor.

Do not use your equipment without a ground connection.

Do not saturate the chamber with material.

Voltage variations can damage electronic components.

If the Plug Cable is damaged a qualified person is needed to change it to avoid damage.

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



(6)

CONFIGURATING A TRIPLE CYCLE.


Within the triple cycle option have the ability to configure.

Three operating temperatures


Three cycle times

Alarm bypass with range of 5 to 30 C.

With the Increment and Decrement keys you can navigate within the (7) menu. To enable or disable each parameter is necessary to select and press the Set button, then the parameter will shade (8) this indicates that we can change its value, with the Increase or Decrease buttons. When we get the desired value of the parameter must be pressed the Set button to add the data in the system. Selecting the option Ejecutar and press Set the cycle started.


(7)


(8)

CANCEL A CYCLE.

To cancel a cycle being executed you have to press the Ok button, followed by the key Cancel, the system ask if you want to cancel the cycle (9) If you pres the Cancel button the cycle is cancel (10) If not pressed after a few seconds the system will eliminate the question and takes us back to the previous state of execution.

(9)

(10)

The care you take in reading and following this instruction will probably determinate the satisfactory service you will receive from your Oven.

UNPACKING.

Carefully remove the Oven from the shipping case. Preserve all paper work for future reference. If damage has occurred from shipment a claim must be filed with the carrier immediately; preserve the shipping container for inspection by the carrier. Contact your dealer or Felisa.

ANALOG MODELS.

FELISA Ovens are available in models with analog control or digital control, both also available with mechanical or gravity convection. The exterior construction is rugged, one/piece cabinet. All models have an operating range of temperature from ambient plus 5 to 250°C. Analog models have a knob to control temperature, digital models are equipped with a digital microprocessor-based temperature control which operates with a pT-100 sensor. It has been designed with high-quality and precision elements, it's totally solid state. All models are equipped with LCD display, with a resolution of 1°C.

INSTALLATION.

Install your Oven using a connection of the adequate voltage. See the label specification before connecting. The Oven can be installed on any surface or table sufficiently firm and strong, leaving a minimum space of 15 cm. Between the Oven and any other vertical surface Keep the area around free of any material to allow ventilation in the area of the bottom. For correct operation it is necessary that the user becomes familiarized with all the available controls and specifications shown in each model.

LOAD IN THE CHAMBER.

It is very important not to overload the chamber. For the right heating the load, it should occupy no more than 3/4 of any dimension of the chamber. The Oven heats up any load by radiation and heat as fast as it absorbs this energy. Try to center the load, the Oven can not be forced without risk of damage.

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You already have a one year warranty for your equipment. If after using this, answer this questionnaire and send it to us within 90 days by fax or email will automatically give one more year warranty!

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CALIBRATION MENU.


The control was calibrated in the plant using certified and reliable patrons, however transportation, temperature or specific location may affect this calibration.

If you have a reliable thermometer, you can calibrate your equipment there against this reference as follows.


Program your oven to the required operating temperature. Place your thermometer in the thermometer holder located on top of the unit, ensuring that the sensor is inside of the chamber. Wait for the temperature to stabilize within the chamber (Approximately 30 minutes after reaching the set temperature. Compare the temperature display against the thermometer and adjust if there are differences as follows.

Return to the start menu (3 and press the Cancel and Increase buttons at the same time, the screens appear (11) and (12) Within this menu we select Escala, and press Set, with the Increment or Decrement keys adjust the temperature in the display and press Ok

If you select the option C. Original, the system will return to the plant calibration values.

(11)

(12)

GRAPHIC AND PRINTING.

From the start menu press the buttons Cancel and Decrement at the same time, you enter to the menu of Graphic and Printing of data (13) In this menu are registered data from the last executed cycle. To observe the generated graph option select Graficar and press OK, by doing this will show the graph on the screen (14) we can return to the previous state by pressing Cancel.

The display shows the option Datos, this refers to the number of data stored in the system memory. To print the stored data, select the option Print and press Ok, then this data is sent to the thermal printer (If your equipment has it

ANALOG OVEN OPERATIONS.

For its operation it is necessary to have on hand a thermometer with a minimum range of 250°C. Placed in the thermometer holder (29), turn on your Oven with the pilot switch (31). Adjust the control knob (33) to the desired temperature. The indicator led (32) will turn on and remain on until reaching the programmed temperature, afterwards it will turn on when the control applies energy to the heating element.

You have to wait 30 minutes approximate in order that the Oven became stable. Verify if the temperature in the thermometer corresponds to that you programmed, if not, readjust the control knob. The experience in operation using the control will allow the operator to identify the exact position for different temperatures. If you are going to work with the same temperature we suggest that once you have the control position for that temperature, turn off the Oven, without moving the control knob and you will always have the same temperature.

DIGITAL OVEN OPERATION.

The control has implemented a Fuzzy algorithm to control the temperature, has a graphic display and a serial communication port for sending data to either a terminal or a thermal printer (optional depending on each model).


BUTTONS DESCRIPTION.

Cancel: Returns you to the previous menu and abort a running cycle.

Set: Selects between the operations of each screen.

Increment: Increases the selected parameter.

Decrement: Decreases the selected parameter.

MENUS.

- 1) Start
- 2) Cycle
- 3) Cycle type
- 4) Configuration
- 5) Print / Graphic
- 6) Calibration.

SPEED CONTROLLER

This device has the power to regulate the RPM's of the motor, offering with this the regulation of air inside the camera. Being 1 its minimum value and 10 maximum of RPM's granted by the engine.

CONTROL ON.

To turn the control on, press the switch to on position. Immediately following will print the screen (1) for a few seconds, the operating system will recognize the operation of the temperature sensor, the strength and integrity of the line voltage. Upon completion of these tasks will hear a beep and the review screen print (2). The control system requires connect to a line grounded to avoid electrical noise problems, also the control system makes detection to line and verify that the supply voltage is correct


(1)


(2)

CYCLE SETTING.

The system can store two different cycles, cycle A, and cycle B. In the cycle menu (3), is where we can make a selection from cycle to use, with the button Set entered the cycle type menu (4), in which we have the option of a single or a triple cycle, with the increment or decrement keys select the desired option and press Set.

(3)

(4)

CONFIGURATION OF A SIMPLE CYCLE.

With the simple cycle you have the option to configure:

Operating temperature

Timer on or off

Cycle time

Bypass alarm with range of 5°C to 30°C.

With the Increment and Decrement keys you can navigate within the (5) menu. To enable or disable each parameter, is necessary to select and press the Set button, then the parameter will shade (6) this indicates that we can change its value, with the Increase or Decrease buttons. When we get the desired value of the parameter must be pressed the Set button to add the data in the system, selecting the option Ejecutar and press Set, the cycle started.