

320 L/min • FE-1800, FE-1801

Presión Garantizada / Guaranteed Pressure	500 mm Hg / 20 in Hg
Desplazamiento en vacío / Free Air Displacement	320.0 L x min / 11.2 cfm
Motor / Motor	1,725.0 rpm
Transmisión / Drive	1 BV / VB
Motor / Motor	120 V
Peso Neto / Net weight	36.0 Kg
Peso Embarque / Shipping Weight	50.0 Kg
Dimensiones W x D x H / Dimensions W x D x H	50.0 x 30.0 x 35.0 cm
Reguladores / Relief valves	Opcional / Optional
Manómetro / Pressure gauge	Opcional / Optional
Vacuómetro / Vacuum gauge	Opcional / Optional
Filtro silenciador / Muffler filter	Estándar / Standard
Cubre banda / Belt guard	Estándar / Standard



FE-1800, FE- 1801
FE-1800L, FE- 1801L

Modelo Model	Desplazamiento // Displacement Litros/min // Liters/min	Vacío máximo Maximum vacuum	Lubricación Lubrication	Velocidad Speed
	mm Hg : 0 : 150 : 250 : 400 : 500	Continua Continuous	Intermitente Intermittent	
FE-1800	: 320 : 250 : 190 : 120 : 60	500 mm Hg	500 mm Hg	No lubricada : 1,440.0 rpm
FE-1800 L	: 320 : 260 : 200 : 130 : 70	500 mm Hg	600 mm Hg	Lubricada : 1,440.0 rpm
FE-1801	: 400 : 320 : 240 : 170 : 90	500 mm Hg	500 mm Hg	No lubricada : 1,725.0 rpm
FE-1801 L	: 400 : 320 : 240 : 180 : 100	500 mm Hg	600 mm Hg	Lubricada : 1,750.0 rpm

700 L/min • FE-1900, FE-1901



FE-1900, FE- 1901
FE-1900L, FE- 1901L

Presión Garantizada / Guaranteed Pressure	500 mm Hg / 20 in Hg
Desplazamiento en vacío / Free Air Displacement	700.0 L x min / 25.5 cfm
Motor / Motor	1,725.0 rpm
Transmisión / Drive	1 BV / VB
Motor / Motor	220 V
Motor / Motor	1½ hp
Peso Neto / Net weight	61.0 Kg
Peso Embarque / Shipping Weight	72.0 Kg
Dimensiones W x D x H / Dimensions W x D x H	60.0 x 40.0 x 40.0 cm
Reguladores / Relief valves	Opcional / Optional
Manómetro / Pressure gauge	Opcional / Optional
Vacuómetro / Vacuum gauge	Opcional / Optional
Filtro silenciador / Muffler filter	Estándar / Standard
Cubre banda / Belt guard	Estándar / Standard

Modelo Model	Desplazamiento // Displacement Litros/min // Liters/min	Vacío máximo Maximum vacuum	Lubricación Lubrication	Velocidad Speed
	mm Hg : 0 : 150 : 250 : 400 : 500	Continua Continuous	Intermitente Intermittent	
FE-1900	: 700 : 550 : 420 : 280 : 140	500 mm Hg	500 mm Hg	No lubricada : 1,060.0 rpm
FE-1900 L	: 700 : 550 : 420 : 280 : 140	500 mm Hg	600 mm Hg	Lubricada : 1,060.0 rpm
FE-1901	: 850 : 700 : 550 : 390 : 230	500 mm Hg	500 mm Hg	No lubricada : 1,325.0 rpm
FE-1901 L	: 850 : 700 : 550 : 390 : 230	500 mm Hg	600 mm Hg	Lubricada : 1,325.0 rpm