Product Specifications



800-253-4942 • IBISCI.com

An ISO 9001 Certified Company

HEPES Free Acid

<u>Product Number</u>	<u>Size</u>
IB01130	50gm
IB01131	250gm
IB01132	500gm
IB01133	1kg

DHVCI	CAL	CDECI	FICATI	ONC
гптэг	CAL	SECT	FICALL	CINO

Appearance	White Crystalline Powder
Molecular Formula	
Molecular Weight	
Ph (0.5M Solution in H ₂ O)	
Purity	
Water	
Heavy Metals	Max. 0.001%
Iron	
Sulfates	Max. 0.05%
DNase assay	None Detected
RNase assay	None Detected
Protease assay	None Detected
Cell Culture Testing	Pass

STORAGE

Store at room temperature. Protect from moisture.

WARNING

Classification: No known hazards

CAS NO.	
CAS#	7365-45-9

RECOMMENDED USE

IBI HEPES, Free Acid is a zwitterionic buffer ideally suited for cell culture in a range of 10 - 25mM. HEPES has a buffering range from 6.8 - 8.2. It is also a component of Tris-HEPES Buffer for polyacrylamide electrophoresis of proteins.