




FE-371 Modelo Analógico




MODELOS ANALÓGICOS / ANALOG MODELS

Termobaños **Felisa®** construidos totalmente en acero inoxidable, su tina de una sola pieza evita fugas y corrosión.

FE-371
60 tubos / 60 tubesFE-372
120 tubos / 120 tubesFE-373
240 tubos / 240 tubes

Features

- Microcontrolled digital control with LCD screen, to observe programmed operating temperatures, cycle time and process status.
- Easy programming in the temperature control now incorporated with "Autotuning" system.
- System with security alarms with high and low temperature limits.
- Fault detection system.
- Timer up to 75 hours.
- Up to 16 temperature profiles and ramp control.
- USB digital communication, through graphic and data acquisition software.
- Analog model with control knob (thermometer required).
- Stainless steel lid and racks for 60 standard size tubes (priced separately).

FE-386
Gradilla 60 tubos / Rack 60 tubesFE-384
Tapa grande / Large coverFE-383
Tapa mediana / Medium coverFE-382
Tapa chica / Small cover